



# The Triple Play

Your Ultimate Hands-Free Decontamination  
& Disinfection System



# RDS 3110

Total Asset Decontamination

# RDS 3110

## ADP-Ex

Designed for Type II Ambulances

ADP-EX Kit includes port, exterior cover and nozzle.



## ADP-PT

Designed for Type I & III Ambulances

ADP-PT Kit includes hose and nozzle assembly.



## ADP-AS

The solutions above require the use of one of our dual-headed hoses.

ADP-AS 8 foot dual-headed hose.



## ADP-RM

Designed for rooms and facilities.

ADP-RM Kit includes port, exterior cover and nozzle.

## Benefits

- Decontaminate even hard-to-reach surfaces
- Consistent and Reliable delivery of disinfectant
- Improves asset turn-around time - back in service faster
- Faster decontamination & disinfection process requiring less labor
- Low cost of operation and ownership
- Safe for sensitive electronic equipment



## Kills Harmful Bacteria, Viruses, and Fungi

Proven effective against:

- Bacteria: MRSA, Staphylococcus aureus, E. coli, Salmonella, Pseudomonas, Legionella
- Viruses: HIV-1, Hepatitis B & C, SARS-CoV-2, Influenza A/Brazil (H1N1), Rotavirus, Norovirus, Poliovirus, Rhinovirus, Human Coronavirus
- Fungi & Mold: Trichophyton mentagrophytes, Aspergillus niger
- Others: Listeria monocytogenes, Klebsiella pneumoniae



## EPA-Approved Disinfectant: Vital Oxide

Active ingredients:

Chlorine Dioxide.....0.200%

Alkyl (60% C14, 30% C16, 5% C12, 5% C18)

dimethyl benzyl ammonium chloride.....0.125%

Alkyl (68% C12, 32% C14) dimethyl ethyl benzyl ammonium chloride.....0.125%

OTHER INGREDIENTS.....99.550%

Total.....100.000%

## APA

### HAND APPLICATOR

Ideal for:

- Exterior settings
- Individual pieces of equipment
- Quick spot-treatment on high contact surfaces



AeroClave, LLC  
4001 Forsyth Road  
Winter Park, FL 32792

Phone : 407-788-3300  
Website : [www.aeroclave.com](http://www.aeroclave.com)  
Email : [info@aeroclave.com](mailto:info@aeroclave.com)

## Decon Large Spaces Effortlessly

### Room Fogging

The RDS 6110 uses an EPA-approved, hospital-grade disinfectant to completely fog large-sized treatment areas such as laboratories, surgical suites, gymnasiums, locker rooms and garages.

### Hand Apply

Use an optional AeroClave Portable Applicator (APA) to apply disinfectant by hand-spraying directly on high contact areas to spot treat patrol cars, ambulances, equipment and exterior spaces.

### Benefits

- Reduces your staffs exposure to dangerous pathogens
- Decontaminate and Disinfect even hard-to-reach surfaces
- Consistent and reliable delivery of disinfectant, not achievable through manual cleaning methods
- Fast turn around and hands-free technology allows you to redeploy staff to highest and best use
- Low cost of operation and ownership
- Safe for sensitive equipment

### Kills Harmful Bacteria, Viruses, and Fungi"

Proven effective against:

- Bacteria: MRSA, Staphylococcus aureus, E. coli, Salmonella, Pseudomonas, Legionella
- Viruses: HIV-1, Hepatitis B & C, SARS-CoV-2, Influenza A/Brazil (H1N1), Rotavirus, Norovirus, Poliovirus, Rhinovirus, Human Coronavirus, Ebola, Mpox, Bird Flu (H5N1)
- Fungi & Mold: Trichophyton mentagrophytes, Aspergillus niger
- Others: Listeria monocytogenes, Klebsiella pneumoniae

### Advanced User Interface

Our sophisticated touchscreen interface allows for minimal training to operate the RDS 6110. Facility and asset management tracking software is included, which will save run data and room configurations. Reports are available in csv file format by download to a portable USB drive.

### Multiple Applications

- EMS/Fire/Police vehicles
- Law enforcement offices, jails, and holding cells
- Hospital facilities & equipment
- Surgical suites, in and out-patient facilities
- Laboratories & research facilities
- \* Nursing homes
- "Assisted Living Facilities
- Veterinary hospitals and animal shelters
- Athletic facilities & sports equipment
- Schools & childcare facilities
- Cruise ships and hotels
- And More



## OPTIONAL ACCESSORIES AND CONSUMABLES



### ADP-PT

Allows for hand-free decontamination and disinfection of the patient compartment of a fire rescue or ambulance.

- Mounting: Installed inside exterior compartment
- Connection: Using ADP connector on RDS 6110
- Hose length: 8 ft



### APA

used for manual hand-application or remote-head fogging.

- Connection: Using ADP connector on RDS 6110
- Hose length: 25 ft
- Tripod: Optional



### Barcode Capabilities

Rooms and assets may be tagged with standard barcodes for easy retrieval of configuration data by the onboard asset tracking software.



# ADS

# AeroClave®



Operate in-station or on the road between transports for maximum flexibility.



Low cost of ownership and operation - about \$1/day.



Finished within 20 minutes to get you back in service.



Environmentally friendly and safe for your staff and equipment.

The modular design and small size allows it to be tucked away in a compartment or hung on the wall




## HIGHLIGHTS

- Fully self-contained system that can operate anywhere
- Fogging provides consistent delivery of disinfectant to all exposed areas, including cabinets and drawers
- I/O for full vehicle integration with on-board control systems (control, monitoring and safety systems)
- Safe for your sensitive electronic equipment
- Remote control for hands-free operation
- All run data available in csv format
- Low cost of ownership and daily operation
- Complete runtime is only 20 minutes
- OEM and after-market installation

## BENEFITS

- Reduces staff and patient exposure to dangerous pathogens
- Consistent and reliable delivery of disinfectant, not achievable through manual cleaning methods
- On-board unit is always ready to be run - On-site or on the go
- Improves unit turn-around time - back in service faster
- Faster decontamination & Disinfection process requiring less labor
- Improved general cleanliness of the vehicle and firehouse
- Decontaminate & Disinfect between transports for maximum protection



Remote Control Included

ADS mounted on cabinet shelf



ADS mounted on compartment wall



# HYDRA SYSTEM



- The Hydra System is the next evolution of AeroClave's patented ADP port technology.
- The Hydra System provides the ability to decontaminate multiple areas or rooms from a single point of operation.
- The entire system is operated by AeroClave's portable RDS 6110 unit, meaning the installed Hydra System requires no power.
- This also provides a portable system that can be used to decontaminate non-Hydra System installed rooms/areas.
- The Hydra System is modular and expandable, allowing for:
  - Multiple Hydra Systems to be installed throughout a facility, if needed.
  - Additional rooms/areas can be added to existing Hydra System Hubs.
  - Additional Hydra Systems can be added in the future.



The Hydra System Hub is the central control point and can connect to up to 12 different area/rooms.

With a low-profile design and no power needed. The Hub can easily be integrated into future or existing facilities.

Two quick-connect adapters supply power, fluid, and air to the system.



The RDS 6110 unit powers and controls the Hydra System Hub, selecting areas for decontamination with precise programmability.



The Hydra System Hub feeds disinfectant fluid and air to remote spray nozzles, supporting up to 12 nozzles that can be installed within a 50-foot range.



## Hydra Applications

- K9 kennels
- Mat Rooms
- Locker rooms
- Weight rooms
- Jail and holding cells
- Administrative offices
- Food service and preparation areas
- Classroom & Patient Room

Spray nozzles feature a tamper-proof, anti-ligature flush design, allowing flexible installation almost anywhere while maintaining safety.






# THE AEROCLAVE PROCESS



“ Our mission is to fight infectious diseases with effective, affordable, and eco-friendly solutions.



### Room Fogging

AeroClave disinfects hard-to-reach surfaces effectively.



### Eco-Friendly

ClO<sub>2</sub> is safer than bleach, with no harmful by-products.



### Destroys Pathogens

Effective against COVID, MRSA, HIV, Ebola, and more.



### Vital Oxide

Broad-spectrum antimicrobial with chlorine dioxide for fogging.



### AeroClave vs. Cleaning

Manual cleaning spreads germs; AeroClave reduces contamination.

### CLO2 DISINFECTION: FAST, EFFECTIVE, AND SAFE

- ✓ ClO<sub>2</sub> stops protein formation, destroying bacteria.
- ✓ Reacts with peptone to starve viruses.
- ✓ Maintains 98% strength for a year, instantly killing microbes.
- ✓ Prevents bacterial mutations and biofilm.
- ✓ Non-irritating, non-corrosive, only eye protection needed.

## PHASES OF THE AEROCLAVE PROCESS

### INJECTION

Aerosolize solution for 5 to 30 minutes.

### DWELL

Allow 10 minutes to distribute

### AERATION

Solution is removed with ventilation methods

Visit our website to download more information on Vital Oxide.

AeroClave, LLC  
4001 Forsyth Road  
Winter Park, FL 32792

Phone : 407-788-3300  
Website : [www.aeroclave.com](http://www.aeroclave.com)  
Email : [info@aeroclave.com](mailto:info@aeroclave.com)

Copyright © 2024 AeroClave, LLC. All rights reserved.






# CASE STUDY

## HOW A DEPARTMENT CUT UNSCHEDULED LEAVE BY 34%

### SICK & UNSCHEDULED LEAVE HOURS COMPARISON



SCAN TO DOWNLOAD CASE STUDY



REDUCING THE INVISIBLE THREAT OF INFECTION IN FIRE & EMS PERSONNEL

AeroClave, LLC  
4001 Forsyth Road  
Winter Park, FL 32792

Phone : 407-788-3300  
Website : [www.aeroclave.com](http://www.aeroclave.com)  
Email : [info@aeroclave.com](mailto:info@aeroclave.com)

## Independent Testing Summary

Disinfection	Contact Time
<b>BACTERIA</b>	
Pseudomonas aeruginosa (ATCC 15442)	10 min
Acinetobacter baumannii (ATCC 19606)	10 min
Staphylococcus aureus MRSA (ATCC 33592)	10 min
Listeria monocytogenes (ATCC 15313)	10 min
Legionella pneumophila (ATCC 33153)	10 min
Salmonella enterica (ATCC 10708)	10 min
Staphylococcus aureus (ATCC 6538)	10 min
Klebsiella pneumoniae (NDM-1) (ATCC BAA-2146)	10 min
Escherichia Coli (ATCC 11229)	10 min
Bordetella bronchiseptica (ATCC 10580)	10 min
Mycobacterium bovis- BCG	10 min
Penicillin- Resistant Streptococcus pneumoniae (ATCC 700677)	10 min
Chlamydia psittaci, Strain 6BC, (ATCCVR-125)	10 min
Streptococcus pyogenes (ATCC 14289)	10 min
Proteus mirabilis (ATCC 7002)	10 min
<b>Disinfection</b>	
<b>VIRUS</b>	
Rotavirus, Strain WA**	5 min
Hepatitis C Virus Bovine Viral Diarrhea Virus**	5 min
Hepatitis B Virus Duck Hepatitis**	5 min
Norovirus Feline Calicivirus**	5 min
Murine Norovirus (MNV-1)**	5 min
Swine Influenza (H1N1)**	5 min
Respiratory Syncytial virus (ATCC VR-26)**	5 min
Human Immunodeficiency Virus (HIV Type 1) Strain HTLV-IIIb**	5 min

Influenza B Virus (Strain B/Taiwan/2/62), (ATCC VR-1735)**	5 min
Hepatitis A Virus (ATCC VR-1541)**	5 min
Hantavirus (Prospect Hill Virus) University of Ontario**	5 min
Canine Distemper Virus, Strain Snyder Hill (ATCC VR)**	5 min
Severe Acute Respiratory Syndrome- Related Coronavirus 2 (SARS-CoV-2) (COVID-19 Virus)**	5 min
Adenovirus 1, Strain Adenoid 71 (ATCC VR-1)**	10 min
Canine Parvovirus, Strain Cornell-780916 (ATCC VR-2016)**	10 min
Canine adenovirus 1, Strain Utrecht (ATCC VR-293)**	10 min
Feline Herpesvirus 1 Strain C-27 (ATCC VR-636)**	10 min
Feline coronavirus, Strain WSU 79-1683 (ATCC VR 989)**	10 min
Canine coronavirus Strain 1-71 (ATCC VR-809)**	10 min
Canine Influenza (H3N8) virus Strain A/Ca/NY/105913/08, Cornell University**	10 min
Feline Infectious Peritonitis Virus, Strain: WSU 79 1146, (ATCC VR-990)**	10 min
Feline panleukopenia virus, Strain Philips-Roxane (ATCC VR-648)**	10 min
Rabies Virus, Strain: Flury (HEP) (ATCC VR-139)**	10 min
Avian Influenza (H5N1) Virus (VNH5N1-PR8/CDC-RG) CDC #2006719965, Charles River Laboratories**	10 min
Herpes Simplex Virus 1 (Human), Strain HF (ATCC VR-260)**	10 min
Porcine Epidemic Diarrhea Virus, Strain CO 2013, USDA APHIS 025-PDV001**	10 min
Porcine Reproductive and Respiratory Syndrome type 2 Virus, Strain NVSL (North America) NVSL 130 PDV**	10 min
Canine herpesvirus, D004 strain (ATCC VR-552)**	10 min
Canine parainfluenza virus, 78-238 strain (ATCC VR-1573)**	10 min
Herpes simplex virus 2, G strain (ATCC VR-734)**	10 min
Human coronavirus, 229E strain (ATCC VR-740)**	10 min
<b>Sanitization</b>	<b>Contact Time</b>
<b>Hard Surface Sanitization No Rinse Required</b>	
Escherichia coli (ATCC 11229)	60 sec
Staphylococcus aureus (ATCC 6538)	60 sec

<b>Sanitization</b>	<b>Contact Time</b>
<b>Non-Food Contact Surface Sanitization No Rinse Required</b>	
Escherichia coli (ATCC 11229)	5 min
Staphylococcus aureus (ATCC 6538)	5 min
<b>Sanitization</b>	<b>Contact Time</b>
<b>To Sanitize Soft Surfaces</b>	
Staphylococcus aureus (ATCC 6538)	5 min
Klebsiella aerogenes (ATCC 13048)	5 min
<b>Sanitization</b>	<b>Contact Time</b>
<b>Carpet Sanitization</b>	
Klebsiella aerogenes (ATCC 13048)	60 min
Staphylococcus aureus (ATCC 6538)	60 min
<b>Fungicide</b>	<b>Contact Time</b>
Trichophyton rubrum (ATCC MYA)	10 min
Trichophyton interdigitale (ATCC 9533)	10 min
<b>Mildewstat</b>	<b>Contact Time</b>
Aspergillus niger (ATCC 6275)	10 min
<b>Insects:</b>	<b>Contact Time</b>
Argentine Ants (Linepithema humile)	10 Min
<b>Organisms not approved by the EPA</b>	
<b>Microorganism</b>	<b>Contact Time</b>
Stachybotrus chartarum	10 Min
Aspergillus fumigatus	10 Min
Alternaria alternata	10 Min
Penicillum sp	10 Min

<b>Mycotoxins:</b>		<b>Contact Time</b>
Aflatoxin		10 Min
Trichothecene		10 Min
Glutoxin		10 Min
Ochratoxin		10 Min
<b>Insects:</b>		<b>Contact Time</b>
Bed bugs (Cimex lectularius)		10 Min
<b>Parasites:</b>		<b>Contact Time</b>
Myxobolus cerebralis (Whirling Disease)		10 Min

Ver 17  
Rev Date: 01/03/2024

#### Section 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product Name Vital Oxide RTU

1.2 SDS No. VO091223

1.3 EPA Reg. No. 82972-1

1.3 Manufacturer/Distributor Vital Solutions, LLC  
3755 Fiscal Ct. Ste. B  
West Palm Beach, FL 33404

1.4 Phone Numbers: (561) 848-1717

1.5 Medical Emergency: (800) 222-1222

Aqueous Oxidant  
HMIS Rating

H	0	Blue
F	0	Red
R	0	Yellow
PE	--	White

#### Section 2: Hazards Identification

**2.1 GHS Classification:** This material is not considered hazardous by the OSHA Hazard Communication Standard 2012 (29 CFR 1910.1200)

**2.2 HMIS Rating:** Health: 0 Flammability: 0 Reactivity: 0 PPE: None

**2.3 Remarks** For full text of Hazard- and EU Hazard-statements: see SECTION 16.

**2.4 Emergency Overview:** Colorless liquid with mild fresh odor. Avoid contact with eyes. Keep out of reach of children.

##### 2.5 Potential Health Effects

- **Eye Contact:** Eye contact may cause mild eye irritation with discomfort.
- **Skin Contact:** Does NOT cause skin irritation and the product is NOT skin sensitizer.
- **Inhalation:** Does NOT cause any respiratory irritation. If consumer product accidentally contacts strong acids in restricted ventilation area, avoid breathing the vapors and allow adequate time for the vapors to disperse before re-entering the restricted area.
- **Ingestion:** Non-Toxic **Carcinogenicity Information:** IARC, NTP, OSHA, and ACGIH list none of the components present in this material at concentrations equal to or greater than 0.1% as carcinogens.

**2.6 Hazard statements:** Eye contact may cause mild eye irritation with discomfort.

##### 2.7 Precautionary statements

- **prevention:** P234 Keep only in original container, or properly labeled second container
- **Precautionary statements – response:** P390 Absorb spillage to prevent waste.

#### Section 3: Composition / Information on Ingredients

Ingredients	CAS Number	Wt %
Oxychlorine Compounds	Mixture	0.200
n-Alkyl Dimethyl Benzyl Ammonium Chloride	68391-01-5	0.125
n-Alkyl Dimethyl Ethylbenzyl Ammonium Chloride	68956-79-6	0.125
Inert Ingredients	Mixture	99.55

At these concentrations none of the ingredients are known to pose any hazards to human health.

## **Section 4: First Aid Measures**

- 4.1 Inhalation:** Does NOT cause any respiratory irritation. If consumer product accidentally contacts strong acids in restricted ventilation area, avoid breathing the vapors, and allow adequate time for the vapors to disperse before re-entering the restricted area.
- 4.2 Skin Contact:** Does NOT cause skin irritation.
- 4.3 Eye Contact:** In case of contact, flush eyes with plenty of water.
- 4.4 Ingestion:** Non-toxic. Give a glass of water.

## **Section 5: Fire Fighting Measures**

- 5.1 Flammable Properties: Flash Point:** Not Available (Non Flammable)
- 5.2 Flammable Limits: Lower Flammable Limit:** Not Established
- 5.3 Burn Rate:** Unknown
- 5.4 Upper Flammable Limit:** Not Established
- 5.5 Flammability Classification:** Non-Flammable liquid
- 5.6 Auto ignition Temperature:** Not Established
- 5.7 Hazardous Combustion Products:** Thermal or other decomposition may yield chlorine dioxide or chlorine.
- 5.8 Extinguishing Media:** N/A (Non-Flammable liquid)
- 5.9 Additional Considerations:** None
- 5.10 FIRE FIGHTING INSTRUCTIONS:** Non-Flammable liquid
- 5.11 NFPA Rating:** Health: 0 Flammability: 0 Reactivity: 0 PPE: None

## **Section 6: Accidental Release Measures**

- 6.1 Spill Clean Up**  
No special cleanup measures are required for the consumer product. To avoid the possibility of "bleaching" the spill should be absorbed with paper towels, and the area rinsed with clean water.
- 6.2 Accidental Release Measures**  
Spills are slippery and should be cleaned up promptly.

## **Section 7: Handling and Storage**

- 7.1 Handling:** Keep away from heat and strong acids. Do not ingest. Keep container closed. Use only with adequate ventilation.
- 7.2 Storage:** Keep container tightly closed and sealed until ready for use. Keep container in a well-ventilated place. Do not store above 90°F (40°C) or near fire or open flame. Store large quantities in buildings to comply with OSHA 1910.106. Do not transfer contents to bottles or other unlabeled containers. Do not reuse empty containers. Keep out of reach of children.
- 7.3 Incompatible materials:** Strong acids
- 7.4 Special Packaging Materials:** None

## **Section 8: Exposure Control/ Personal Protection**

- 8.1 Engineering Controls:** Use in adequately ventilated areas.
- 8.2 Personal Protective Equipment:**
- 8.3 Eye/Face Protection:** Not required for consumer product.
- 8.4 Skin Protection:** Not required for consumer product.
- 8.5 Respirators:** None required for normal use. If consumer product accidentally contacts strong acids in restricted ventilation area, Avoid breathing the vapors, and allow adequate time for the vapors to disperse before re-entering the restricted area.
- 8.6 Exposure Limits:**
- **Oxychlorine Compounds**
  - **n-Alkyl Dimethyl Ethylbenzyl Ammonium Chloride:**
    - PEL (OSHA): Not available PEL (OSHA): Not available
    - TLV (ACGIH): Not available TLV (ACGIH): Not available

**• n-Alkyl Dimethyl Benzyl Ammonium Chloride:**

PEL (OSHA): Not available

TLV (ACGIH): Not available

**Section 9: Physical and Chemical Properties**

<b>Appearance:</b>	Colorless liquid	<b>Odor:</b>	Mild-Fresh
<b>Physical State:</b>	Liquid	<b>pH:</b>	8 - 9
<b>Boiling Point (°F):</b>	212	<b>Solubility in Water:</b>	100%
<b>Freezing Point (°F):</b>	32	<b>Vapor Pressure (mm Hg):</b>	Not Available
<b>Volatile Organic Compounds (VOC):</b>	None	<b>Evaporation Rate:</b>	Less than Ether
<b>Specific Gravity:</b>	1.003 @ 68°F (20°C)	<b>Density (lb/gal):</b>	8.40 @ 68°F (20°C)

**Section 10: Stability and Reactivity****10.1 Chemical Stability:** The product is stable.**10.2 Incompatibility with other Materials:** Strong acids**10.3 Conditions to avoid:** Contact with strong acids**10.4 Hazardous Polymerization:** Will not occur.**10.5 Hazardous Decomposition Products:** Thermal or other decomposition may yield chlorine dioxide or chlorine.**Section 11: Toxicological Information****11.1 TOXICITY TESTING**

- **Acute Inhalation** – Studies with Wistar Albino rats exposed to a respirable aerosol made from a solution of Vital Oxide at a level of 2.08 mg/l for four hours resulted in no deaths and no abnormal necropsy observations.
- **Eye Contact** – Studies with New Zealand white rabbits showed this product is a very mild ocular irritant; mild conjunctival irritation was observed, but cleared within 24 hours.
- **Skin Contact** – Study of dermal toxicity in New Zealand white rabbits showed the product to be non-toxic: Dermal LD<sub>50</sub>> 5,000 mg/kg of body weight; Study of dermal irritation in New Zealand white rabbits showed the product is not a dermal irritant. In Dermal Sensitization studies, Vital Oxide was determined not to be a sensitizer.
- **Swallowing** - Acute oral toxicity in albino rats: Non-toxic LD<sub>50</sub>>5,000 mg/kg of body weight.

**11.2 EPA TOXICITY RATING – IV** This is the lowest category on the scale and is designed for substances that are the least hazardous.**Section 12: Ecological Information****12.1 Environmental Hazards:** No data available.**12.2 Environmental Fate:** No data available.**Section 13: Disposal Considerations****13.1 Waste Disposal:** Treatment, storage, transportation, and disposal must be in accordance with applicable Federal, State/Provincial and Local regulations.**Section 14: Transport Information****14.1 Shipping Information:** Not regulated by DOT, IMO/IMDG and IATA/ICAO for ground, air or ocean shipments.

---

---

## Section 15: Regulatory Information

---

### 15.1 U.S. Federal Regulations:

- **TSCA:** All components appear in TSCA Inventory
- **OSHA:** Refer to Section 8 for exposure limits.

### 15.2 CERCLA SARA Hazard Category:

- **Section 311 and 312:** This product has been reviewed according to the EPA "Hazard Categories" promulgated under Sections 311 and 312 of the Superfund Amendment and Reauthorization Act of 1986 (SARA Title III) and is considered, under applicable definitions, to meet the following categories: None
- **Section 313:** This product contains following substances subject to the reporting requirements of Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 and 40 CFR Part 372: None.

**15.3 State regulations:** State Right to Know information is not provided.

**15.4 California prop. 65 (no significant risk level):** None

**15.5 International Regulations:**

**15.6 Canadian WHMIS:** Not controlled.

**15.7 Canadian Environmental Protection Act (CEPA):** Additional information available upon request.

**15.8 EU Regulations:** Additional information available upon request.

**15.9 EPA/FIFRA: Registration number: 82972-1**

This chemical is a pesticide product registered by the United States Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS), and for workplace labels of non-pesticide chemicals. The hazard information required on the pesticide label is reproduced below. The pesticide label also includes other important information, including directions for use.

---

---

## Section 16: Other Information

---

The information is furnished without warranty, expressed, or implied, except that it is accurate to the best knowledge of Vital Solutions, LLC. The Data on this sheet related only to the specific material designed herein. Vital Solutions, LLC, assumes no legal responsibility for the use or reliance on this data.

<b>Issue Date:</b>	12 September 2023
<b>Revision Date:</b>	0
<b>Revision Note:</b>	0

---

---

End of SDS