

Automobile Decontamination

ADP-Ex



The ADP-Ex is designed for those vehicles without an exterior compartment, but need a cost-effective option to be used in conjunction with AeroClave RDS units to decontaminate using AeroClave's hands-free disinfectant system.



Benefits

- Reduces your staff's exposure to dangerous pathogens
- Fast turnaround and hands-free technology allows you to redeploy staff to highest and best use
- Consistent and reliable delivery of disinfectant, not achievable through manual cleaning methods
- Extremely easy to use with minimal training required
- Fully utilize the capabilities of your RDS unit
- Safe for sensitive electronic equipment
- Low cost of operation; approximately \$1/run

Components

The ADP-Ex is installed on the outside of the vehicle and a nozzle assembly is mounted in the interior patient compartment. These are linked with two lengths of color coded tubing that is easily fed through the structure of the vehicle.

ADP-RM SPECIFICATIONS

- Maximum treatment area: 700 cubic ft
- RDS compatibility: RDS 3110 and RDS 6110
- Kit includes:
 - ADP Panel Mount
 - Dimensions: (Dia) 3.0in x (D) 2.5in
 - Material: Polypropylene
 - Connection: AeroClave's ADP plug-and play
 - Exterior Cover
 - Dimensions: (L) 4.75in x (W) 3.75in x (D) 1.25in
 - Type: Weather-Resistant
 - Material: Thermoplastic Glass-Filled
 - Nozzle Assembly
 - Dimensions: (L) 1.5in x (W) 1.5in x (D) 2.25in
 - Material: Stainless steel and plastic
- Two 15 ft 1/4 inch polyurethane lengths
- Warranty: 12 months parts and labor
- Training: Available online or on-site



Overview

The ADP-Ex is designed for those vehicles without an exterior compartment. Utilize your RDS unit to fog the interior of the ambulance, while eliminating excess manual labor and exposure in contaminated areas. AeroClave's automated process ensures consistent operation and results, disinfecting even those hard-to-reach surfaces.



Required Hose

This solution requires the use of one of our quad-headed hoses.

Part Number: ADP-Ex

How does it work?

Connect one end of the ADP-AS hose to the RDS and the other to the ADP-Ex on the vehicle. Open all interior cabinets and close the exterior doors. Set the RDS to AMBULANCE MODE and push the start button. The process virtually eliminates the human element of disinfecting, allowing for consistent results. Best of all, it's all done in less than 20 minutes.

ADP-Ex CONNECTION



Compatible with RDS systems

Every AeroClave RDS unit has at least one ADP port as standard equipment. Optional models are equipped with up to three (3) ADPs for even more flexibility.

Installation

Easily installed in under 30 minutes. The kit comes with everything you need to install as an OEM, after-market or DIY installer.



AeroClave, LLC
4001 Forsyth Road
Winter Park, FL 32792
Phone: 407-788-3300
Web: www.aeroclave.com

Decontamination Technology